

Renishaw CMM Products Division

PRODUCT BULLETIN – PBC1780

Product:	MSP3	Status:	Open
Title:	MSP3 obsolescence notice	Date:	11th March 2015

Distribution:	Internal	Subsidiaries	OEMs	Distributors	Retrofitters
----------------------	----------	--------------	------	--------------	--------------

Originator: Gareth Puntis	
Reviewer: Andy Holding	

Summary:

This bulletin is to give notice that the MSP3 manual scanning probe will be made obsolete from March 2015. Renishaw has been making this product since January 2008.

The MSP3 manual scanning probe was developed for the Quantum FaroArm made by FARO Technologies. The product overcome the need to press the capture buttons on the arm, ensuring that data was only taken when the stylus was in contact with the part.

It is Renishaw's policy to support an obsolete product, through repair, RBE, or replacement for five years. This policy will be applied to the MSP3 starting April 2015.

Part numbers:

The affected part numbers are as follows:

Part number	Description
A-3053-2131	MSP3 probe
A-3053-3200	Renishaw MSP3 kit
A-3053-2200	FARO MSP3 kit
A-3053-2300	MSP3 NDI probe

If you have any questions or comments on this bulletin please visit www.renishaw.com/cmmsupport and use the "Contact" facility. This will ensure that your call is logged and processed efficiently.